FLEX
Monoclonal Mouse
Anti-Human
Nucleophosmin
Clone 376
Ready-to-Use
(Dako Autostainer/Autostainer Plus)

Code IS652

ENGLISH	
Intended use	For in vitro diagnostic use. FLEX Monoclonal Mouse Anti-Human Nucleophosmin, Clone 376, Ready-to-Use, (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. This antibody is useful for the identification of both wild-type nucleophosmin (NPM) and the mutated counterpart NPMc+. Nucleophosmin is localized in the nuclei of normal cells in multiple tissues, but is aberrantly accumulated in the cytoplasm of leukemic blasts in a large subgroup of acute myeloid leukemia cases with a normal karyotype (1-3). The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonyms for antigen	B23, NO38, NPM, numatrin (1).
Summary and explanation	Nucleophosmin (NPM) is a phosphoprotein belonging to the protein family of nucleoplasmins (1). NPM is predominantly localized in the nucleus of cells in many tissues (4). However, NPM has the ability to shuttle between the nucleus and cytoplasm, and this shuttling activity and its proper subcellular localization might be crucial for cellular homeostasis (1). The NPM gene (<i>NPM1</i>) is frequently overexpressed, mutated, deleted or rearranged in many human malignancies, and, depending on its expression level and gene dosage, NPM seems to function as either an oncogene or a tumor suppressor. Complete functional loss of NPM leads to aneuploidy and activation of cell proliferation checkpoints (1). Mutations in exon 12 of the <i>NPM1</i> gene is causal for generation of a mutant NPM counterpart, NPMc, which is aberrantly localized in the cytoplasm of leukaemic blasts cell in a large subgroup of primary acute myeloid leukemia (AML) cases with a normal karyotype (1-4). NPMc display a multilineage distribution in AML specimens of all subtypes except M3 (acute promyelocytic leukemia), M4eo (acute myelomonocytic leukemia with eosinophilia) and M7 (acute megakaryocytic leukemia) (2, 4). It is likely that mutation and cytoplasmic accumulation of NPMc, called NPMc+, is a primary leukemogenic event, as it may interfere with the normal NPM function. NPMc is not observed in any primary AML cases with a specific genetic abnormality, in secondary AML specimens or hematopoietic or extrahematopoietic neoplasms other than AML (2, 4). NPMc is only observed in primary AML with a normal karyotype (2, 4, 5). Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.
Reagent provided	Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L sodium azide. <u>Clone</u> : 376. <u>Isotype</u> : IgG1, lambda.
Immunogen	Recombinant full length NPM protein.
Specificity	In Western blotting of U937 cell lysates, AML cell line lysates (2), and lysed OCI-AML3 cells (3), the antibody labels a band of 37 kDa corresponding to NPM.
Precautions	1. For professional users. 2. This product contains sodium azide (NaN ₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. 3. As with any product derived from biological sources, proper handling procedures should be used. 4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. 5. Unused solution should be disposed of according to local, State and Federal regulations.
Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation including materials required but not supplied	The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. Pre-treatment with heat-induced epitope retrieval (HIER) is required using Dako PT Link (Code PT100/PT101). For details, please refer to the PT Link User Guide. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). <u>Paraffin-embedded sections</u> : Pre-treatment of formalin-fixed, paraffin-embedded tissue sections is recommended using the 3-in-1 specimen preparation procedure for Dako PT Link. Follow the pre-treatment procedure outlined in the package insert for EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Code K8010/K8004). Note: After staining the sections must be dehydrated, cleared and mounted using permanent mounting medium. <u>Deparaffinized sections</u> : Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link and following the same procedure as described for paraffin-embedded sections. After staining the slides should be mounted using aqueous or permanent mounting medium. The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended.
Staining procedure including materials required but not supplied	The recommended visualization system is EnVision™ FLEX, High pH, (Dako Autostainer/Autostainer Plus) (Code K8010). The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols: Template protocol: FLEXRTU2 (200 µL dispense volume) or FLEXRTU3 (300 µL dispense volume) Autoprogram: NPM (without counterstaining) or NPMH (with counterstaining) The Auxiliary step should be set to "rinse buffer" in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to "none". This ascertains comparable wash times. All incubation steps should be performed at room temperature. For details, please refer to the Operator's Manual for the dedicated instrument. If the protocols are not available on the used Dako Autostainer instrument, please contact Dako Technical Services. Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. If the evaluating pathologist should desire a different staining intensity, a Dako Application Specialist/Technical Service Specialist can be contacted for information on re-programming of the protocol. Verify that the performance of the adjusted protocol is still valid by evaluating that the staining pattern is identical to the staining pattern described in "Performance characteristics".

(117051-003)IS652/EFG/MNI/2009.12.04 p. 1/4

Dako Denmark A/S | Produktionsvej 42 | DK-2600 Glostrup | Denmark | Tel. +45 44 85 95 00 | Fax +45 44 85 95 95 | CVR No. 33 21 13 17

	Counterstaining in hematoxylin is recommended using EnVision™ FLEX Hematoxylin, (Dako Autostainer/Autostainer Plus) (Code K8018). Non-aqueous, permanent mounting medium is recommended. Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include liver and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics" in all positive specimens. The recommended negative control reagent is FLEX Negative Control, Mouse, (Dako Autostainer/Autostainer Plus) (Code K8018) (Code IS750).
Staining interpretation	Cells labeled by the antibody display either nuclear or cytoplasmic and nuclear staining.
Performance characteristics	<u>Normal tissues</u> : The antibody labels wild-type NPM in nuclei of gut, intestine, kidney, liver, lung, lymph node, pancreas, prostate, spleen, thymus and tonsil. Nuclear labeling is observed in all cells in cerebellum but large neurons display additional staining in the cytoplasm. The hepatocytes in the liver show a moderate to strong nuclear staining reaction, whereas the scattered mitotic and post-mitotic cells show a weak to moderate cytoplasmic and nuclear staining reaction. <u>Abnormal tissues</u> : The antibody labeled both the wild-type NPM and mutant NPMc proteins (2, 3) and resulted in both nuclear and cytoplasmic labeling of leukemic blast cells in 208/591 (~35%) primary AML specimens, and in 142/230 (~61%) of primary AML specimens with a normal karyotype. In contrast, nuclear and cytoplasmic labeling of leukemic blast cells was only observed in 24/263 (~9%) AML cases with an abnormal karyotype, and none of these cases represented specific genetic abnormalities commonly observed in primary AML. Nuclear labeling alone was observed in 135/135 (100%) secondary AML specimens and in 980/980 hematopoietic or extrahematopoietic neoplasms other than AML (2).

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